

We bring more variety to the Vitrex® Sterilance Press II safety lancet product line



**Vitrex® Sterilance Press II
30G Sterile Safety lancet**

GENTLE: Maximum patient comfort

Fine gauge, tri-bevel tip for virtually painless sampling

Single drop blood volume

Colour Code	Configuration	Penetration Depth	Article No.	Application
Green	● 18G blade	1.8 mm	393213	Special baby care and maternity units for heel puncture
Yellow	● 26G needle	1.8 mm	394213	Diabetes clinics/hospital wards for blood glucose monitoring
Orange	● 21G needle	2.2 mm	395213	Haemoglobin (blood banks), cholesterol, blood typing, coagulation, blood gas and other diagnostic tests
Pink	● 21G needle	2.8 mm	396213	As Orange but for thicker adult skin
Blue	● 21G needle	1.8 mm	397213	As Yellow but for increased blood flow
Purple	● 28G needle	1.8 mm	398213	Low blood flow
Turquoise	● 30 G needle	1.8 mm	399213	Low blood flow – less pain NEW ITEM

Pre-Order NOW

Questions about this product?
Call us 039201 702055 or email info@labdiasys.de

Labor Diagnostik Systeme GmbH

Am Obstgarten 10
D-39326 Wolmirstedt

Phone +49 39201 70 20 55
Fax +49 39201 70 96 99

info@labdiasys.de
www.labdiasys.de

Testing become simple